# Z7 - Studding Magnetic Relaxation of Protons by Spin Echo Method

Physics Laboratory II

Faculty of Physics, Astronomy and Applied Computer Science, Jagiellonian University

The goal of the exercise is an investigation of pulse methods of nuclear magnetic resonance (NMR). Student learns how NMR spectroscope works, plans and makes measurements (also prepares chemical solutions) and this way gets competences in laboratory work. The measuring of relaxation times  $T_1$  and  $T_2$  of protons in water solutions of organic and inorganic compounds is a basic result of exercise. Additionally there are also measurements of temperature and concentration dependence of  $T_1$  and  $T_2$ .

### Preparatory questions

- 1. NMR classical description [1]
  - Magnetic moment in magnetic field equation of motion in laboratory and in rotating frame
  - Effective field
  - Bloch equation
- 2. NMR quantum description [1]
- 3. Pulse methods for investigation of NMR [1]
  - Resonance pulses of radio frequency  $\pi$  and  $\pi/2$
  - Free Induction Decay
  - Spin echo
  - Different pulse sequences
- 4. Nuclear magnetic relaxation phenomenon [1]
  - Spin-net relaxation  $T_1$
  - Spin–spin relaxation  $T_2$
- 5. Investigation of temperature dependence activation energy [1]

## **Computational assignments**

Calculate magnetic nuclear resonance of hydrogen, deuterium and helium in a field of induction  $2\mathrm{T}$ 

### Apparatus and materials

The set-up (Figure 1) consists of:

- Pulse spectrometer PS15 (2) controlled by PC (1)
- Electromagnet (7) with electrical supplier (5) and cooling system
- Nuclear stabilizer of magnetic field JTM83 (6)
- Heating-cooling system (4, 8) with a temperature controller (3)

### Experiment

- Principle of operation and maintenance of PS15 spectrometer with external electromagnet creating B<sub>0</sub>, stabilized by JTM83
- Measuring of  $T_1$  and  $T_2$  for protons in vicinity of Cu ions in water solution of  $CuSO_4$  in different temperatures using appropriate pulse technique
- Measuring of  $T_1$  and  $T_2$  for protons in water solutions of proteins in different temperatures using appropriate pulse technique (especially for biophysics)
- Measuring of  $T_1$  and  $T_2$  for protons in distilled water in different temperatures using appropriate pulse technique

#### Exploration of the set-up:

- 1. Switching on the water cooling of electromagnet. It is recommended to wait a quarter of an hour.
- 2. Switching on: PC, electromagnet electrical supplier, spectrometer PS15, JTM83 stabilizer
- 3. Launching SPI program (lock-off function)
- 4. Changing a current of electromagnet to achieve a resonance value of B0 (FID signal on a display)
- 5. Activation of stabilization of magnetic field
- 6. Finding time duration of pulses  $\pi$  and  $\pi/2$
- 7. Applying appropriate pulse sequences establish  $T_1$ and  $T_2$  for a sample in a room temperature and below (with a help of electrical cooling system)

#### Data analysis

The experimental data should fit the appropriate theoretical dependencies:

- Finding  $T_1$  and  $T_2$  of protons in CuSO4 in different temperatures
- Finding  $T_1$  and  $T_2$  of protons in distilled water solution in different temperatures
- Finding activation energy for water





Figure 1: Apparatus

# Safety rules

Be careful with strong magnetic field and liquid nitrogen. Student with peacemakers should avoid the exercise. Keep watches, phones, memory sticks etc. away from magnetic field.

### References

[1] J.P.Hornak, The Basics of NMR

